

A pregnant woman with long brown hair, wearing a black short-sleeved dress and black high-heeled boots, is standing by a large window. She is leaning forward, tending to a tall, dark vase filled with purple flowers. The room is bright and modern, with a wooden chair and a large wooden vase on a rug in the foreground. The window looks out onto a green landscape.

medi

Venous Health catalog

Compression solutions



Who is medi?

medi® is committed to helping people all around the world live a more independent, productive and satisfying life while managing circulatory, lymphatic and ambulatory issues. medi takes its position as a global leader in medical compression and orthopedics seriously by investing in research, education and innovation by providing the latest technologies and enabling our patients to not only manage their challenges, but to enjoy life to its fullest.

“I feel better”

is a statement of our commitment to our customers, clinicians and patients to deliver upon our promises of unparalleled service, educational support and healthy outcomes in the clinical communities we serve.

The medi brand collection

mediven®

duomed®

CIRCAID®



CEP

MESI

protect

restiffic®

REJUVA®

Advancis®



4 compression therapy

6 mediven technologies

- 12 mediven sheer & soft
- 13 mediven comfort
- 14 mediven comfort sculpt
- 15 mediven comfort vitality
- 16 mediven plus
- 17 mediven forte
- 18 mediven angio
- 20 mediven for men classic & select
- 21 mediven active
- 22 mediven dual layer
- 23 mediven size charts

24 duomed

- 25 duomed transparent
- 26 duomed advantage
- 27 duomed relax & patriot

28 Rejuva

30 CEP

32 circaid

34 medi compliance solutions

Keys to successful compression therapy

compression & compliance

To ensure best patient outcomes, medi products are made with expansive sizing options to ensure a proper fit and delivery of measurable compression levels.

Quality and Standards

Manufactured to exacting standards, medi compression products rise above the competition in rigorous testing to provide accurate and effective medical compression.

Listed Compression	Product	Size	Tested Ankle Compression
15-20 mmHg	sheer & soft	III	18.3 mmHg
15-20 mmHg	comfort	III	18.4 mmHg
15-20 mmHg	for men	III	17.9 mmHg
20-30 mmHg	sheer & soft	I	24.4 mmHg
20-30 mmHg	comfort	IV	24.9 mmHg
20-30 mmHg	for men	II	24.2 mmHg
20-30 mmHg	plus	VI	25.3 mmHg
30-40 mmHg	sheer & soft	II	33.3 mmHg
30-40 mmHg	comfort	V	35.1 mmHg
30-40 mmHg	for men	VI	34.6 mmHg
30-40 mmHg	plus	IV	35.1 mmHg
30-40 mmHg	forte	VII	36.3 mmHg

- FDA registered
- ISO certified medical device manufacturer

Clinically Proven

In a published clinical study, 70% of patients wearing mediven stockings were still wearing them after 16 months.

Accurate compression dosage

- Perfect fit sizing for accuracy across the broad spectrum of patients
- Adjustable compression wraps with Built-in-Pressure System™ (BPS™)

Comfort through quality

- Acu-Knit technology for soft, smooth compression
- Advanced fibers that provide moisture regulation
- Anti-microbial properties to control odor

Easy application and care

- Soft elastic knit in foot
- Juxtaposed bands for alternative garment applications
- Machine wash and dry



ready-to-wear product lines



custom mediven



circaid

comprehensive & confidence

medi's individualized solutions fit every lifestyle with our expansive product options. Our comprehensive product lines are always at the cutting edge of development:

ready-to-wear mediven product lines:

- fashion - sheer & soft, comfort, comfort vitality, comfort sculpt, for men classic and select, active and harmony
- medical - plus, forte, dual layer, angio and thrombexin

custom mediven:

- flat knit - 550, cosy 450, mondi 350
- circular knit - comfort, plus, forte, angio

circaid:

inelastic circaid wraps - juxtalite® & juxtalite® hd, juxtafit®, juxtacures®, graduate™ and reduction kit™, circaid® profile

lifestyle compression products:

- Rejuva - fashionable compression legwear
- CEP - compression sportswear

PCS:

- brio & genius

medi compression quality

30-day wear guarantee

The patient can return the product for any reason related to the fit, comfort or style of the medi product.

6-month compression guarantee

medi compression products will deliver the stated compression for a minimum of 6 months.

customer service

- With a 98% fulfillment rate, receive timely and accurate processing of your orders.
- medi's certified, knowledgeable customer service and sales staff are available to answer all of your questions.



mediven[®]

innovative compression for life

Compression therapy can be challenging for those with therapeutic and lifestyle needs. The comprehensive offering of mediven products uses the latest technologies in materials and knitting design to construct compression products that help overcome fit and application challenges.

Perfect Fit technology

Should these women be fit with the same size stockings?

Precision has been the foundation of mediven stocking since their inception. The secret is our ability to consistently make a wide range of sizes with specific dimensions to create a nearly customized fit.

IF THEY WERE:

5'2"
Julie



— 8 3/8"
actual size ankle
circumference

Would both women
get a good fit?

Would both women
receive accurate
compression at the ankle?

Would Lisa find her
stockings hard to put on?

Could Julie find
her stockings falling
down and bunching
at the ankles?

5'8"
Lisa



9 3/4"
actual size ankle
circumference

The answer is of course, **NO!**

These women should not be fit with the same size stockings. That's why medi offers the most complete range of sizes available today.

With medi ready-made stockings, all your patients receive the proper, custom-like fit and precise compression that the doctor prescribed, everytime.

mediven


other brands

ankle sizes	7	5
length options	2	sometimes
ankle size range	3/4"	1 1/4" and up



Acu-Knit technology

The knitting process in mediven products creates a unique open mesh pattern called Acu-Knit. This delivers a multitude of therapeutic and compliance benefits to enhance the patient experience with the goal of creating lifelong customers.



competitor (left) vs mediven plus (right)
at 40x magnification



Why choose mediven Acu-Knit?

Compression integrity = maintains compression throughout the day and the lifetime of the garment

Compression profile = smooth and accurate compression keeps the recommended gradient compression level

Comfort = Acu-Knit allows air to recirculate, creating a natural feel of freshness all day

Durability = Acu-Knit creates stronger materials to reduce tears and punctures, without compromising therapy



Clima Fresh

Clima Fresh


The integrated Clima Fresh function has an antibacterial action to prevent odor formation. Bacteria are kept away from the knit and the fiber-deep hygiene is maintained.



without Clima Fresh

Clima Comfort

Moisture is transported away from the skin to allow for greater temperature control and breathability.



fabric

skin surface

Additional technologies found in mediven products:



Anatomical Fit The Anatomical Fit feature found in many of our unisex and men products allow for greater foot size variance and provide enhanced comfort.



mediGrip Top Band Elastic stockings have two-way stretch that allows for comfort during mobility. The mediGrip top band helps secure calf-length stockings in place without slippage.



Soft Elastic Foot Since not all venous health patients require edema control in the foot area, some mediven products contain a soft elastic foot with reduced compression made for easier application and long-term comfort.

Easy to care and maintenance makes for greater compliance : Most mediven products can be machine washed and dried. Although this may not be a technology, it is a BIG contributing factor when patients make purchasing decisions and influences their wearing frequency!!



- Hand wash or machine wash only using gentle detergents on delicate cycle at 100°F.
- Do not use fabric softeners, chlorine based or Woolite™ type detergents.
- Most mediven models are suitable for machine dryers using the gentle cycle. Do not dry over heaters, in direct sunlight or by wringing out the garments.
- Do not dry clean or iron your medi garments
- Use the medi stocking guardian wash bag to wash your compression garments with extra care and to prolong their lifespan.
- Follow the specific wear and care instructions provided with each product.





calf stocking

mediven sheer & soft

Timeless stocking for the fashionable woman.

Indications:


- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

Features:

- medi's most transparent compression for fashion-conscious patients.
- Clima Comfort ensures unparalleled softness, enhanced breathability and moisture regulation.
- Clima Fresh keeps stockings odor-free and fresh.
- Open and closed toe designs provide a reinforced foot for superior durability.
- Acu-Knit textile technology delivers all day therapy for ease of application and a custom-fit feeling.



panty



open toe detail

attractive lace top band



mediven sheer & soft

styles				
calf	thigh w/lace silicone top band	panty	maternity panty	
colors available				
WHEAT	NATURAL	TOFFEE	CHARCOAL	EBONY
				NAVY
Not all colors available in every style-please refer to dealer price list.				
compression levels available: 8-15, 15-20, 20-30, 30-40 mmHg				



Also available in petite



Also available in open toe

mediven comfort


The most comfortable stocking in its class.

Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

Features:

- Clima Comfort ensures unparalleled softness, enhanced breathability and moisture regulation.
- Clima Fresh keeps stockings odor-free and fresh.
- Semi-transparent fabric helps conceal blemishes with confidence.
- Soft, luxurious feel provides all-day comfort year round.
- Acu-Knit textile technology delivers all day therapy for ease of application and a custom-fit feeling.



comfort panty top



double-knitted comfort sole

mediven comfort

styles	
calf	extra-wide calf
thigh	thigh with lace or beaded silicone top band
panty	maternity panty

colors available

WHEAT	SANDSTONE	NATURAL	CHOCOLATE	NAVY	EBONY
-------	-----------	---------	-----------	------	-------

Not all colors available in every style-please refer to mediusa.com.

compression levels available: 15–20, 20–30, 30–40 mmHg

silicone top band with beaded and lace options



Also available in petite

Also available in open toe





medical compression legging

mediven comfort sculpt

Designed for comfort; built to shape.

Indications:


- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

Features:

- Opaque material in fashionable, classic colors.
- Tummy shaper sculpts waist and supports lower back.
- Flexible fabric with buttery soft feel.
- Added support improves circulation and reduces leg pain and discomfort.
- Soft, luxurious fibers wick moisture away.
- Full coverage in ankle for comfort and quality compression.



tummy shaper




flexible fabric

mediven comfort sculpt

styles



colors available



compression levels available: 15-20 mmHg

Also available in petite



full coverage

mediven comfort vitality

The most comfortable stocking in its class - with pattern.

Indications:

- Minor to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

Features:

- Fashionable, Herringbone-pattern for men and women.
- mediGrip technology keeps sock in place, thanks to a seamless silicone top band.
- Clima Fresh keeps stockings odor-free and fresh.
- Clima Comfort blended with Lanolin keeps legs comfortably dry all day.
- Available in standard and petite lengths.
- Acu-Knit textile technology delivers all day therapy for ease of application and a custom-fit feeling.



mediven comfort vitality

styles	colors available
 <p>calf</p>	<div style="display: flex; justify-content: space-around;"> <div style="background-color: #4b3621; color: white; padding: 5px;">CHOCOLATE</div> <div style="background-color: #1a1a1a; color: white; padding: 5px;">EBONY</div> <div style="background-color: #333; color: white; padding: 5px;">CHARCOAL</div> </div> <p>compression levels available: 15–20, 20–30, 30–40 mmHg</p>

 Also available in petite



Herringbone-ribbed pattern



mediGrip

Herringbone-ribbed pattern detail





thigh-length stocking with knitted border




men's leotard



Thigh with waist attachment

silicone top band



mediven plus

Broad selection with enhanced containment.

Indications:

- Moderate to severe varicose veins.
- Heavy, achy legs.
- Healed ulcers.
- Early-stage lymphedema.
- Moderate to severe edema.
- Ideal for patients with skin issues (ulcers and staph infections).

Features:

- Opaque coverage conceals blemishes and varicose veins with confidence.
- Strong, durable material for superior support.
- Clima Comfort provides unmatched softness and breathability, keeping legs fresh.
- Clima Fresh keeps stockings odor-free and fresh.
- Combination of both elastic and inelastic properties for serious edema containment.
- Acu-Knit textile technology delivers all day therapy for ease of application and a custom-fit feeling.



mediven plus

styles					colors available
calf	calf w/beaded silicone top band	extra-wide calf w/beaded silicone top band	thigh w/beaded silicone top band	panty	<div style="background-color: #d9c8a1; padding: 5px; margin-bottom: 5px;">BEIGE</div> <div style="background-color: #000000; color: white; padding: 5px; margin-bottom: 5px;">BLACK</div> <p>Not all colors available in every style—please refer to mediusa.com.</p>
<i>only in open toe</i>					
extra-wide calf	thigh	thigh with waist attachment	maternity panty	men's leotard	



Also available in petite



Also available in open toe

mediven forte

Comfortable, secure stocking for edema control.

Indications:



- Moderate to severe varicose veins.
- Heavy, achy legs.
- Moderate to severe edema.
- Healed ulcers.
- Early stage lymphedema.
- Helps prevent deep vein thrombosis (DVT).

Features:

- Opaque coverage conceals blemishes and varicose veins with confidence.
- Our strongest material ensures maximum all-day support.
- Clima Comfort provides exceptional breathability and temperature control.
- Clima Fresh keeps stockings odor-free and fresh.
- Combination of both elastic & inelastic properties for serious edema containment.
- Strongest elastic mediven material available.
- Open toe only.



mediven forte

styles			
			
calf	extra-wide calf	extra-wide calf with beaded silicone top band	thigh with beaded silicone top band
colors available			
CARMEL			
compression levels available: 30-40, 40-50 mmHg			

 Also available in petite



calf stocking with open toe



thigh length stocking with open toe



open toe detail



calf stocking

mediven angio

With
Merino
plush!

Safety in treatment - down to the toes.

Indications:

- Chronic venous disease.
- Chronic venous disease and mild to moderate PAD.
- Chronic venous disease and diabetes mellitus Type 2/diabetic PNP.

Features:

- Merino plush in the entire foot provides unique padding throughout foot and ankle to reduce pressure on bony prominences.
- Extra-long extended flat seam reduces pressure points on toes and feet.
- Solid rib structure provides the ideal ratio between resting and working pressure.
- Knitted top band features a broad design with a secure hold.



merino plush (shown inside out)



extra long flat-seam (shown inside out)

mediven angio

styles


calf

colors available

CARAMEL BLACK

compression levels available: 15-20, 20-30 mmHg

Also available in petite



solid ribbed structure



Do you see these patients?

- 200MM people globally have PAD¹
- 422MM people globally have Diabetes²
- Approximately 40% of the population has chronic Venous insufficiency


Guide to compression therapy for patients with peripheral arterial disease or diabetes.



Study proves mediven[®] angio is safe compression for mild-moderate PAD

Key findings

- mediven angio is safe for CVI patients with mixed PAD or mixed Diabetes Mellitus
- Microcirculation measured stable in all subjects using 15-20 and 20-30 mmHg
- Improved perfusion was seen showing the benefit of compression stockings





calf sock



classic pattern showing anatomically shaped foot with right and left markings

select pattern showing the extra-wide top band

mediven for men classic & select

For dress up or dress down days.

Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

Classic Features:

- Classic dress sock design for a professional look.
- Lanolin blend provides premium softness and feel.
- Anatomical left and right foot design for the perfect fit.
- Clima Fresh keeps stockings odor-free and fresh.
- Acu-Knit textile technology delivers all day therapy for ease of application and a custom-fit feeling.

Select Features:

- Casual pattern for work and leisure.
- Soft, plush fabric matches your relaxed lifestyle.
- Soft, extra-wide mediGrip top band for a comfortable, stay-put fit.
- Anatomical left and right foot design for the perfect fit.
- Clima Fresh technology inhibits bacteria growth, keeping legs clean and odor-free.
- Acu-Knit textile technology delivers all day therapy for ease of application and a custom-fit feeling.



only for select pattern

mediven for men classic & select

styles				
calf	calf tall	extra-wide calf	extra-wide calf tall	
colors available				
classic	TAN	GREY	NAVY	BROWN
				BLACK
Not all colors and styles available in the 8-15mmHg offering - please refer to mediusa.com .				
select	TAN	GREY	BLACK	
compression levels available: 8-15, 15-20, 20-30, 30-40 mmHg				

Features mediGrip

T Also available tall

mediven active

Cushioned sock for active wear.

Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).






Features:

- Energize and maintain healthy legs with a sport sock for active lifestyles.
- Cushioned foot delivers a stable and secure fit for freedom of movement.
- Performance-based compression profile provides incredible energy return.
- Clima Comfort blend provides exceptional breathability, keeping legs dry during activity.
- Zinc-rich fibers release odors and inhibit bacteria growth, perfect for all-day freshness.
- Acu-Knit textile technology delivers all day therapy for ease of application a a custom-fit feeling.



Clima Comfort Perfect Fit Variety Acu-Knit Cushioned Sole

mediven active

styles	
 calf standard	 calf tall
colors available	
 WHITE	 GREY
 BLACK	
compression levels available: 15–20, 20–30 mmHg	

T Also available tall



calf sock



articulating flex zone



cushioned sole



closed toe liner and open toe outer stocking



soft top band on liner and outer stocking



All mediven dual layer stocking systems include an easy to use donning aid!

mediven dual layer stocking system

Two-layer system for active wound healing.

Indications:

- Active venous ulcers.
- Healed venous ulcers.
- Prevention of reulceration.

Features:


- dual layer system provides accurate compression for effective healing.
- Easy application
- One standard length fits men & women.
- Under liner is soft and flexible.
- Kit includes 2 liners and 2 outer stockings.
- Designed to work with a primary dressing.
- Clima Fresh keeps stockings odor-free and fresh.
- Ideal for self management of venous ulcer treatment.



Relevant Codes:

- HCPCS A6531 gradient compression stocking, below knee, 30-40 mmHg, each. Report with correct modifiers, RT (right) and/or LT (left).
- HCPCS A6532 gradient compression stocking, below knee, 40-50 mmhg, each. Report with correct modifiers, RT (right) and/or LT (left).
- Code usage should be verified by responsible parties in accordance with all published guidelines.

styles



colors available

BEIGE

liner only

WHITE

compression levels available: 30-40 and 40-50 mmHg

mediven dual layer stocking system 30-40 and 40-50 mmHg* closed toe

*When both liner and stocking are worn together.

mediven size charts

mediven sheer&soft, comfort, vitality, comfort sculpt, active, for men classic & select 15-20, 20-30, 30-40 mmHg size chart

circumference	I	II	III	IV	V	VI	VII
panty & thigh	16½" – 22½" 42–57 cm	17¾" – 24¼" 45–61.5 cm	19¼" – 26¼" 49–66.5 cm	20½" – 28¾" 52–73 cm	22" – 31" 56–79 cm	23¾" – 33¼" 60–84.5 cm	25½" – 34½" 64–87.5 cm
extra-wide calf	N/A	N/A	16½" – 19½" 42–49.5 cm	17¾" – 20¾" 45–53 cm	19" – 22" 48–56 cm	20¼" – 23¼" 51.5–59 cm	21½" – 24½" 54.5–62 cm
calf	11" – 14" 28–35.5 cm	11¾" – 15¼" 30–39 cm	12¾" – 16½" 32.5–42 cm	13¾" – 17¾" 35–45 cm	14½" – 19" 37–48.5 cm	15½" – 20¼" 39.5–51.5 cm	16¼" – 21½" 41.5–54.5 cm
ankle	7" – 7¾" 18–20 cm	8" – 8¾" 20.5–22.5 cm	9" – 9¾" 23–25 cm	10" – 10¾" 25.5–27.5 cm	11" – 11¾" 28–30 cm	12" – 13" 30.5–33 cm	13¼" – 14¼" 33.5–36 cm
leg length	calf ▶ standard: 15"–19" (38–47.5 cm) petite: 12–15" (31–38 cm) panty & thigh ▶ standard: 28"–34" (71–86.5 cm) petite: 23"–28" (58.5–71 cm)						


mediven plus, angio and forte size 20-30, 30-40, 40-50 mmHg size chart

circumference	I	II	III	IV	V	VI	VII
panty & thigh	15" – 22½" 38–57 cm	16" – 24¼" 41–62 cm	17¼" – 26¼" 44–67 cm	18½" – 28¼" 47–72 cm	22" – 31" 56–79 cm	21" – 32¼" 53–82 cm	22¼" – 34¼" 56–87 cm
thigh w/silicone top band	16½" – 19¼" 42–49 cm	17¾" – 20¾" 45–53 cm	19" – 22¼" 48–57 cm	20½" – 24" 52–61 cm	22" – 25¾" 56–65 cm	23½" – 27½" 60–70 cm	25" – 29" 64–74 cm
extra-wide calf w/ silicone top band w/out silicone top band	13¼" – 15" 34–38 cm	14½" – 16½" 37–42 cm	15¾" – 18" 40–46 cm	17" – 19½" 43–50 cm	18" – 21" 46–54 cm	19¼" – 22¾" 49–58 cm	20½" – 24½" 52–62 cm
calf w/ silicone top band w/out silicone top band	10¼" – 13¼" 26–34 cm	11" – 14½" 28–37 cm	11¾" – 15¾" 30–40 cm	12½" – 17" 32–43 cm	13¼" – 18" 34–46 cm	14¼" – 19¼" 36–49 cm	15" – 20½" 38–52 cm
ankle	6½" – 7¼" 16.5–18.5 cm	7½" – 8¼" 19–21 cm	8½" – 9¼" 21.5–23.5 cm	9½" – 10¼" 24–26 cm	10½" – 11¼" 26.5–28.5 cm	11½" – 12½" 29–32 cm	12¾" – 13¾" 32.5–35 cm
leg length	calf ▶ standard: 15" – 19" (38–47.5 cm) petite: 12" – 15" (31–38 cm) panty & thigh ▶ standard: 28" – 34" (71–86.5 cm) petite: 23" – 28" (58.5–71 cm)						

mediven dual layer stocking system size chart

circumference	S	M	L	XL	XXL
calf	11" – 15" 28–38 cm	12¼" – 16½" 31–42 cm	13½" – 18" 34–46 cm	14½" – 19½" 37–50 cm	18" – 22" 46–56 cm
ankle	7" – 8½" 18–21.5 cm	8½" – 9¾" 21.5–25 cm	9¾" – 11" 25–28 cm	11" – 12½" 28–32 cm	12½" – 14½" 32–37 cm
leg length	standard: 15" – 19" (38–47.5 cm)				

mediven sheer & soft 8–15 mmHg thigh and panty



mediven sheer & soft 8–15mmHg calf

select by shoe size

S	4-5
M	5½–7½
L	8–10
XL	10½–12



duomed
value you can trust



duomed by medi

The duomed compression line delivers effective venous disease management in a comprehensive range of colors and styles, all backed by the medi product guarantee. duomed is the full solution provider in the compression therapy segment that checks all of the boxes of quality, service, selection and price.

duomed transparent NEW

A sheer, fashionable compression stocking for women.

Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent Deep Vein Thrombosis (DVT).

Features:


- Guaranteed graduated compression for improved leg health.
- Sheer material for a fashionable look.
- Lightweight material for comfortable all-day wear.
- Thigh garments with comfortable lace silicone top band.
- Fashionable styles to match any outfit.
- Smooth, breathable fiber blend is easy to apply.



duomed transparent calf in nude




duomed transparent thigh in nude

styles		
 calf	 thigh with lace silicone topband	 panty
colors available		
NUDE	BLACK	
compression levels available: 15–20, 20–30 mmHg		

 Also available in petite

transparent 15-20, 20-30 mmHg size chart

Leg Circumference	S	M	L	XL	XXL
 Panty	20-26" (51-66 cm)	22-28" (56-71 cm)	24-30" (61-76 cm)	26-31.5" (66-80 cm)	27.5-32.35" (70-82 cm)
CG Thigh	19.25-24.5" (49-62 cm)	21.25-26.75" (54-68 cm)	23.75-29" (60-74 cm)	25.5-30.25" (65-77 cm)	27.5-32.35" (70-82 cm)
CC Calf	11-15" (28-38 cm)	12-16.5" (31-42 cm)	13.5-18" (34-46 cm)	14.5-19.75" (37-50 cm)	17.5-21.5" (44-54 cm)
CB Ankle	7.5-8.75" (19-22 cm)	8.5-9.75" (22-25 cm)	9.75-11" (25-28 cm)	11-12.25" (28-31 cm)	12-13.5" (31-34 cm)
	Leg Length		Standard		Petite
	calf		15-17.5" (39-44 cm)		13-15" (34-38 cm)
	thigh & panty		28.25-32.75" (72-83 cm)		24.25-28" (62-71 cm)



duomed transparent panty in black



duomed • 25

duomed advantage

A soft, opaque stocking for men and women.

Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent Deep Vein Thrombosis (DVT).

Features:

- Guaranteed graduated compression for improved leg health.
- Opaque, durable material to hide blemishes.
- All open toe garments come with donning aid.
- Thigh garments designed with comfortable beaded silicone top band.










duomed advantage calf in almond



duomed advantage thigh in black



duomed advantage panty in beige

styles			
			
calf	thigh with beaded silicone topband	panty	maternity panty
colors available			
			
compression levels available: 15-20, 20-30, 30-40 mmHg			




Also available in petite



Also available in open toe

advantage 15-20, 20-30, 30-40 mmHg size chart

Leg Circumference	XS	S	M	L	XL	XXL
	17.75" – 23.75" 45–60 cm	20-26" (51-66 cm)	22-28" (56-71 cm)	24-30" (61-76 cm)	26-31.5" (66-80 cm)	27.5-32.35" (70-82 cm)
	N/A	19.25-24.5" (49-62 cm)	21.25-26.75" (54-68 cm)	23.75-29" (60-74 cm)	25.5-30.25" (65-77 cm)	27.5-32.35" (70-82 cm)
	N/A	N/A	N/A	17-20" (43-51 cm)	18-22" (46-56 cm)	20-24" (51-61 cm)
	10.25" – 13.5" 26–34 cm	11-15" (28-38 cm)	12-16.5" (31-42 cm)	13.5-18" (34-46 cm)	14.5-19.75" (37-50 cm)	17.5-21.5" (44-54 cm)
	6.75" – 7.5" 17–19 cm	7.5-8.75" (19-22 cm)	8.5-9.75" (22-25 cm)	9.75-11" (25-28 cm)	11-12.25" (28-31 cm)	12-13.5" (31-34 cm)
	Leg Length				Standard	Petite
	calf				15-17.5" (39-44 cm)	13-15" (34-38 cm)
	thigh & panty				28.25-32.75" (72-83 cm)	24.25-28" (62-71 cm)

duomed relax

Cushioned compression sock for casual footwear.


Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent Deep Vein Thrombosis (DVT).

Features:

- Guaranteed graduated compression for improved leg health.
- Cushioned foot bed ensures all-day comfort.
- Opaque material to hide blemishes.
- Soft, flexible fabric is easy to apply.

styles



calf

colors available

BLACK WHITE

compression levels available:
15-20, 20-30 mmHg




duomed relax in black



duomed relax in white

relax & patriot size chart

Leg Circumference	S	M	L	XL	XXL
	11-15" (28-38 cm)	12-16.5" (31-42 cm)	13.5-18" (34-46 cm)	14.5-19.75" (37-50 cm)	17.5-21.5" (44-54 cm)
	7.5-8.75" (19-22 cm)	8.5-9.75" (22-25 cm)	9.75-11" (25-28 cm)	11-12.25" (28-31 cm)	12-13.5" (31-34 cm)
				Standard	
	Leg Length			15-17.5" (39-44 cm)	

duomed patriot

The classic, ribbed dress sock for men.


Indications:

- Mild to moderate spider veins.
- Slight to severe varicose veins.
- Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent Deep Vein Thrombosis (DVT).

Features:

- Guaranteed graduated compression for improved leg health.
- Opaque material to hide blemishes.
- Classic ribbed design seamlessly matches any outfit.
- Comfortable, breathable material for all-day wear.
- All colors available in all styles & compression levels.

styles



calf

colors available


NAVY BLACK TAN

compression levels available:
15-20, 20-30, 30-40 mmHg



duomed patriot in black

detail of ribbed pattern





REJUVA[®]

Stylish compression



The Rejuva portfolio of products is a fashion-forward line of compression legwear that complements products from the mediven line. These fun medical support socks have a wide variety of designs, prints, colors and patterns that add a pop of flair to any wardrobe while also providing functional therapy.

REJUVA®




NEW Crew



Multiple colors and compression levels available.

15-20 mmHg
20-30 mmHg



SIZING BY LEG CIRCUMFERENCE	S	M	L	XL
A (ANKLE)	7-8.5	8.5-10	10-11.5	11.5-13
B (CALF)	11-15	12-17	13-18	14-20
C (THIGH)	16-24	18-28	21-31	23-32
l (INSEAM)	Up to 29	Up to 31	Up to 32	Up to 33

IN INCHES



EMBRACE THE PRESSURE

Today's market leader in athletic compression legwear was first embraced in the endurance sport industry and has been adopted across many sport genres for performance, recovery and injury prevention. CEP products may be used by some of today's greatest athletes but are highly recommended for every skill level. Regardless of the sport or the ability, CEP has the right product for every athlete.

TALL COMPRESSION SOCKS 4.0



CALF SLEEVES 4.0



MID CUT SOCKS 4.0



LOW CUT SOCKS 4.0



NO SHOW SOCKS 4.0



SIZE CHARTS

CALF CIRCUMFERENCE	IN INCH	9.5-12	12.5-15	15.5-17.5	18-20
	IN CM	25-31	32-38	39-44	45-50
WOMEN		II	III	IV	-
MEN		-	III	IV	V



Measure your calf circumference at the widest point.



CIRCAID[®]

Easy, Measurable, Effective

circaid is a full line of inelastic compression therapy wraps for treatment of venous health and lymphatic conditions. The circaid mascot is the giraffe, the tallest animal in the world. Yet, due to its inelastic skin, it does not suffer from venous or edema conditions.

The juxtalite products make up the venous health segment of the circaid portfolio. The juxtalite products are a perfect alternative to elastic compression due to the juxtaposed application straps. This ease of application is why juxtalite is the most popular compression wrap in wound care.

circaid juxtalite & juxtalite hd lower leg

The proven inelastic compression wrap.

juxtalite Indications:

- Mild edema
- CVI
- Active & healed venous ulcers

juxtalite hd Indications:

- Moderate edema
- CVI
- Active & healed venous ulcers
- Phlebolympheoedema

Features:

- Multiple compression levels (20-30, 30-40, 40-50 mmHg).
- juxtalite band system allows for easy application and instant adjustability without removal.
- Built-In-Pressure system (BPS™) guide card provides targeted compression ranges for prescribed dosage.
- Clima Fresh prevents static, odor and propagation of bacteria in the garment.
- Perfect fit with 8 sizes & 2 lengths.
- 6-month warranty.
- Machine washable.

Includes:


- circaid juxtalite and juxtalite hd system includes 1 legging, 1 pair of compressive undersocks, and 1 Built-In-Pressure System (BPS) guide card
- circaid power added compression band™ (pac band™) or juxtalite ankle foot wrap can be purchased separately.



circaid juxtalite with pac band




circaid juxtalite hd with compressive undersocks



circaid juxtalite & juxtalite hd

circumference in centimeters								
size	S	M	M - X full calf	L	L - X full calf	XL	XL - X full calf	XXL
cC	26 - 36	33 - 44	44 - 54	41 - 51	51 - 61	37 - 48	48 - 61	53 - 64
cB	19 - 22	22 - 27	22 - 27	25 - 32	25 - 32	30 - 37	30 - 37	35 - 42
actual product length in centimeters								
size	short				long			
∠B - D	28				33			
leg length		short			long			
∠ floor - D1		<44			>44			


Product lengths may differ from body measurements.





medi compliance solutions

At medi, compression therapy is a comprehensive solution. Our medi compliance solutions assist our patients with application, skin care and the laundering of garments to create the most positive patient experiences using our products.



medi butler
item # 90002

medi butler

Features:

- Optimal 20 degree angle for ease of use.
- Promotes patient compliance and independence.
- Time saving
- Protects garment from misuse



medi "big" butler
Item # 90000

medi "big" butler

- Larger diameter allows for extra-wide calf circumferences.



assure butler
item # 90007

medi assure butler

- The economical donning option.



medi butler off

- The medi butler off enables easy removal of your compression stockings by extending your reach.



easy on foot slipper

- Application aid (comes standard with all mediven open toe products).
- For use with open-toe products and arm sleeves.
- One size fits all.



2 in 1 application aid

- Donning and doffing aid for compression stockings and arm sleeves.
- For open and closed toe garments
- Hard-wearing and tear-resistant
- Light and very compact



medi day gel

Benefits

- Fresh feeling all-day long.
- Reduces the feeling of tension and skin sensitivity.
- Contains horse chestnut for firmer vessel walls.
- Refreshing menthol and witch hazel stabilizes the tissue structure.
- Available in 50 mL tubes.

Use day gel directly before applying compression garments.



medi night creme

Benefits

- For a comforting recovery over night.
- Enhances the regeneration of distressed skin.
- Jojoba oils provide a soothing action on the outer skin layer.
- Available in 50 mL tubes.

Apply night creme after taking compression garments off for a relaxing night's sleep.



it stays!®

- This adhesive securely and comfortably holds knee-high stockings, thigh-high stockings and arm sleeves in place.
- it stays! is odorless, hypo-allergenic and easy to wash away.



medi textile gloves

- Sizes: small, medium, large and X-Large.
- Available in individual pairs and in boxes of 12.



medi care™ wash

- medi care wash is an advanced detergent, specially formulated to clean mediven stockings and sleeves.
- This pure, organic and biodegradable solution is specially formulated to safely remove accumulated body oils and perspiration from garment fibers.



mediven compliance kit

- Contains stocking application gloves, medi care wash and easy care wash bag.



Please join our medi USA
channel on Youtube
medi YouTube Channel
[https://www.youtube.com/user/
mediUSAvideos](https://www.youtube.com/user/mediUSAvideos)



Please join us on Twitter
<https://twitter.com/ifeelbettermedi>
@ifeelbettermedi



Please join us on
[https://www.facebook.com/
mediusaifeelbetter](https://www.facebook.com/mediusaifeelbetter)
@mediusaifeelbetter



Please join us
[https://www.linkedin.com/
company/medi-usa](https://www.linkedin.com/company/medi-usa)



Please join us on [https://www.instagram.com/
mediusaifeelbetter/](https://www.instagram.com/mediusaifeelbetter/)
@mediusaifeelbetter

medi USA
6481 Franz Warner Parkway
Whitsett, NC 27377
T 800-633-6334
F 888-570-4554
mediusa.com